

# Designed for **efficiency** in your lab.

# UniCel DxH 600 Coulter Cellular Analysis System

Beckman Coulter is once again leading the way in lean laboratory practices and efficiency with the introduction of the UniCel DxH 600 Coulter Cellular Analysis System.

Delivering exceptional quality of results, first-pass accuracy and efficiency through cutting-edge technology the UniCel DxH 600 Coulter Cellular Analysis System delivers unparalleled value for your mid-to-large volume laboratory—all with a small footprint.




Blood Banking
Centrifugation
Chemistry
Flow Cytometry
Hematology

Immunoassay Information Systems Lab Automation Molecular Diagnostics Rapid Diagnostics

#### **Superior Quality of Results**

- Multi-Angle Flow Cytometric Digital Morphology technology and advanced algorithms ensure exceptional quality of results
- 29 measurements per cellular event to accurately characterize WBC differential samples, ensuring first-pass accuracy and efficiency
- 2D-dataplots, surface plots and 3D-cube formats allow for enhanced visual cellular investigation




#### **Improved Efficiency**

- Extended Decision Rules capability and Auto-Rerun and Reflex testing reduce TAT and the need for manual operator intervention
- Run standard and micro-volume
   Whole Blood Sample Tubes in
   Automatic Cassette Mode
- Single Aspiration Pathway, simplifies calibration and quality control processes, while eliminating mode-to-mode procedures
   BD Microtainer MAP Microtubes\* utilized in DxH Cassette
- Intelligent Workload Tracking, software allows for workload tracking by day, hour and test

### **Extended Testing Capabilities**

- Multi-Angle Flow Cytometric Digital Morphology technology delivers Automatic NRBC enumeration, without the need for additional reagents or expensive stains.
- Automatic Reticulocyte and Body Fluid Analysis capabilities



**Data Plots and Surface Plots for NRBCs and Reticulocytes** 

#### **Standardization of Process**

- Same technology and User Interface as the UniCel DxH 800 allow for true standardization of process and results across the DxH series and reduces laboratory personnel training needs
- Shared consumables with the DxH 800 allow for simplified inventory management



Easy-to-Use User Interface Reduces Training

## UniCel DxH 600 Coulter Cellular Analysis System Specifications

#### **Available Test Menu**

CBC WBC, UWBC, RBC, HGB, HCT, MCV, MCH, MCHC, RDW, RDW-SD, PLT, MPV

NE, LY, MO, EO, BA, NRBC, NE#, LY#, MO#, EO#, BA#, NRBC#

Retic RET, RET#, MRV, IRF

Body Fluids (Spinal, Serous and Synovial) TNC, RBC

#### **Throughput**

Differential

Up to 100 samples/hr

#### **Aspiration Volumes**

Open/Closed Vial 165 µL

#### **Power**

90-264 VAC and 48-62 HZ

DxH 600 System with System Manager and printer draws 1,185 watts DxH 600 System with System Manager, Printer and UPS draws 2,390 watts

Weight and Dimensions	Depth	Width	Height	Weight
Specimen Processing Module	78.7 cm	76.2 cm	96.5 cm	~117.5 kg
	(31.0 in)	(30.0 in)	(38.0 in)	(260 lb)
Pneumatic Supply	51.0 cm	21.6 cm	36.8 cm	~18.4 kg
	(20.0 in)	(8.5 in)	(14.5 in)	(40.5 lb)

#### Whole Blood Performance — Measuring Range

WBC	$0.050 - 400.000 \text{ X } 10^3 \mu\text{L}$
RBC	0.005 – 8.500 X 10 <sup>6</sup> μL
HGB	0.10 - 25.50 g/dL
MCV	50.00 - 150.00 fL
PLT	$3.0 - 3,000.0 \times 10^{3} \mu$ L
NRBC	0.00 - 600.00/100 WBC
Retic	0.000 - 30.000 %

#### **Data Management**

Computer Microsoft Windows\*\* XP

Database 40,000 Patient results with Graphics


QC with Levy-Jennings Graph, XB/XM for Moving

Averages, Daily Check, Intelligent Quality Monitoring

#### **Sample Transport Capabilities**

Capacity 20 five-tube cassettes
Bar Codes Digital Bar Code with 2D Bar Code capability
Sample ID Up to 22 characters

#### www.beckmancoulter.com/DxH



Intelligent Workload Tracking Software allows for workload tracking per day, hour and even test



<sup>\*</sup>BD and Microtainer MAP Microtubes are trademarks of Becton, Dickenson and Company

<sup>\*\*</sup>Microsoft and Windows are trademarks of Microsoft Corporation

Beckman Coulter, the stylized logo, Coulter and UniCel are trademarks of

Beckman Coulter, Inc. and are registered with the USPTO. DxH is a trademark
of Beckman Coulter, Inc.